



Technical Data Sheet

- 51465 Disposable Bariatric Repositioning Sling, 1200# (1 pkg. - 1 sling)

Device Name	Disposable Bariatric Repositioning Sling
Device Description	Patient repositioning sling for use with a patient transfer lift and a two point spreader bar to reposition or turn a patient in bed. Intended for single patient use. It is not washable.
Unit Description	1 Disposable Bariatric Repositioning Sling
Units per Package	1
Package Dimensions (in)	Approximately 17 L x 11 W x 6 H (1 pkg. - 1 sling)
Package Weight (lbs)	4
Unit Dimensions (in)	84 L x 48 W
Unit Weight (lbs)	4
Lifting Capacity (lbs)	1200
HIBC Code	M66851465
Materials	Body: Breathable Cotton-Polyester blend Strapping: Polypropylene Thread: Polyester
Radio Translucent	Yes
Health/Safety	Contains no latex or PVC
Recyclable	No, combustible waste
Packaging	Minimal – all packaging materials are recyclable
FDA Device Class	Class 1
FDA Regulatory Number	880.5510 (component of)
European Device Class	Class 1
European Regulation	Medical Device Directive 93/42/EEC
GMDN Code	40535



QNET BV
Hommerterweg 286
6436 AM Amstenrade
The Netherlands

