

# TOTALIFT

Lateral Transfer  
Transport Stretcher  
Treatment Chair

WY'EAST  
MEDICAL



The Totalift is 3 devices in 1: a lateral transfer device, a transport stretcher and a treatment chair. The Totalift moves patients laterally in a controlled and safe way from bed to stretcher, no manual lifting required. It can then convert easily from a horizontal stretcher to a comfortable chair.



Early Mobility



Transport



Non-Powered




90600



Patient Tray\*



IV / O<sub>2</sub> Storage



Upper Side Rails\*  
& Headrest\*



Pressure  
Management  
Mattress  
Upgrade\*

\*optional accessories

## Specifications

- Maximum Load 400 pounds
- High Position 36.5"
- Low Position 21"
- Width 27"
- Wheel base 24" x 23.5"



1-18005-R00

## Features

- Non-powered/low maintenance
- Tuck-away lower side rails
- Adjustable leg and foot rest
- Built-in O<sub>2</sub> and IV pole holders
- Comfort Mattress
- Early Mobility tool
- Cardiac, Neuro and ICU applications
- No sling required
- Alternative option to patient lifts

800.255.3126 | [wyeastmedical.com](http://wyeastmedical.com)