

SUMMIT⁵⁰⁰

Powered Sit-to-Stand

WY' EAST
MEDICAL



Early Mobility

The Summit is a moderate-assist, powered sit-to-stand. Its removable foot plate and IV pole/O₂ tank holder make it a must-have for your Early Mobility initiatives. This device has one of the smallest footprints on the market, making maneuverability easy in tight spaces.



Transfers



Edge of Bed



Therapy



51750




Patient Scale*



Removable Footplate



IV Pole/ O₂ Tank Holder*



Leg Straps*

*optional accessories

Features

Small footprint for easy storage and maneuverability in tight spaces

Inter-sling compatibility

O₂ tank/IV pole holder

QR coded user manual

Removable footplate for patient ambulation

Integrated patient scale option

Specifications

- | | |
|--------------------------------------|----------------|
| • Maximum Load | 500 lb |
| • Number of Lifts per Charge | 35-40 |
| • Turning Diameter | 46.1" Diameter |
| • Lifting Speed Unloaded | 1.4 in/s |
| • Leg Under bed/Stretchers Clearance | 2.5" or 3.4" |
| • Standard Door Width Clearance | 23.8" |



1-18002-R01

800.255.3126 | wyeastmedical.com