

# VISTA400

Non-Powered Stand Aid

WY' EAST  
MEDICAL



The Vista is a stand aid designed with the caregiver and patient in mind. This device has color-coded levers and knobs for quick identification and ease of use. Its removable foot plate makes it a tool for safe ambulation as well as a minimal assist transfer device.




**Transfers**



**Maternity**




**Therapy**



**Removable Footplate**



**Pneumatic Stand Assist**



**O<sub>2</sub> Tank/IV Pole Holder\***



**Swing-Away Seat**

\*optional accessories

## Features

- Removable footplate
- Unique, pneumatic sit-to-stand assist
- QR coded user manual
- Non-powered, always ready for use
- Small footprint for easy storage and maneuverability in tight spaces
- Seat pads rotate up and/or out to accommodate varying patient sizes
- Base opens and closes
- O<sub>2</sub> tank/IV pole holder

## Specifications

- Maximum Load 400 lb
- Turning Diameter 46.1" Diameter
- Leg Under bed/Stretchers Clearance 2.5" or 3.4 "
- Standard Door Width Clearance 23.7"



1-18003-R01

800.255.3126 | [wyeastmedical.com](http://wyeastmedical.com)