

# DENALI600

Power-Driven Mobile Lift

WY'EAST  
MEDICAL



Car Extraction

The Denali is a power-driven mobile lift designed with the functionality of a ceiling lift. The ergonomic design of this lift eliminates the push/pull forces, reducing the risk of injury to caregivers. The height adjustable mast and sling bar strap operate independently for easy car extractions and floor recovery.



Floor Recovery



Patient Transfers



Therapy

## Specifications

- |                                       |                  |
|---------------------------------------|------------------|
| • Maximum Load                        | 600 lb           |
| • Number of Lifts per Charge          | 200              |
| • Turning Diameter                    | 53" Diameter     |
| • Lifting Speed Unloaded              | 1.75 in/s        |
| • Leg Under bed/Stretcher Clearance   | 2.5", 3.4", 4.6" |
| • Standard Door Width Clearance       | 28"              |
| • Small or Standard Sling Bar Options | 15" or 24"       |




1-18004-R01




52300



Power-Drive  
Wheels



Intuitive Controls



Mesh Bag with  
Wipeable Sling\*



O<sub>2</sub> Tank/IV Pole  
Holder

\*optional accessories

## Features

- Power driven to eliminate push/pull forces
- Car extraction in any size vehicle
- O<sub>2</sub> tank/IV pole holder
- Removable hand control for patient ambulation/therapy
- Low bed clearance base option
- Ceiling lift alternative
- Base opens and closes
- Sling made of friction-reducing material, designed for car extractions
- Easy floor recovery in tight spaces

800.255.3126 | [wyeastmedical.com](http://wyeastmedical.com)