

PINK SLIP ELITE

Friction-Reducing Device

WY'EAST
MEDICAL



Boosting

The Pink Slip Elite is a friction-reducing device designed with the caregiver in mind. The Pink Slip Elite is single-patient-use and can be used multiple times for many patient handling tasks with the same patient, eliminating cross-contamination. It is specially formulated with a proprietary, durable woven nylon anti-friction treatment.



Pivoting



In-Bed Exercises



Lateral Transfer

Specifications

- Load Recommendations
- Unit Dimensions
- Use Cycle

See reverse side
33" x 57"
Single patient,
multiple use



1-18012-R00

631016: 1 sleeve of 25

631018: Case of 100



Unfolded On Bed

Features

Contains no latex or PVC

Radiolucent

Recyclable

Warm-touch technology

800.255.3126 | wyeastmedical.com

LOAD RECOMMENDATIONS





Caregiver Force Required to Initiate Movement

Liberty Mutual Group loss prevention field organization has been a leader in evaluating lifting, lowering, pushing and carrying tasks. Liberty Mutual has developed metrics for evaluating these tasks with the goal of controlling employee disability costs associated with the manual handling operations. An important conclusion from the published guidelines is: "As a general rule of thumb, designing manual tasks for greater than 75% of the female work population will offer the best protection from manual handling injuries".

Manual Materials Handling Guidelines, 2012, Liberty Mutual Insurance. 55 Lbs initial pulling force once every 30 minutes.

This graph shows the pulling forces required by the caregiver increases progressively for the draw sheet as the patient weight increases. Using the Pink Slip Elite allows the caregiver to safely move 253% more weight when compared to a draw sheet, thus reducing the risk of injury.



This graph shows the pulling forces required by the caregiver increases progressively for the low-density slider as the patient weight increases. Using the Pink Slip Elite allows the caregiver to safely move 60% more weight when compared to using a low-density slider, thus reducing risk of injury with the Pink Slip Elite.