

# DISPOSABLE SLINGS

Loop-Style



## Specifications

## Features

Highback



- **One-Size-Fits-Most**
  - Weight Capacity: 600 lb (272 kg)
  - Product #: **51561**/VA#: **51590**
  - Sold in cases of 5

One size, fits the average-sized patient, improving compliance and inventory management

Built-in head support

Intuitive, single loop straps

Quick application, 4-strap sling

\*S, L, XL sizes available in reusable option

Repositioning



- **Standard**
  - Weight capacity: 1000 lb (454 kg)
  - Dimensions: 42" x 66"
  - Product #: **51411**/VA#: **51591**
  - Sold in cases of 5
- **Long**
  - Weight Capacity: 1000 lb (454 kg)
  - Dimensions: 42" x 84"
  - Product # **52092**/VA#: **52082**
  - Sold in cases of 5
- **Bariatric**
  - Weight Capacity: 1200 lb (544 kg)
  - Dimensions: 48" x 84"
  - Product # **51465**/VA#: **51592**
  - Sold each

3-loop straps for head positioning

2-loop straps on body of sling for sizing and comfort

**Breathability Study:** Independent heat & water vapor, body index study performed on fabric for breathability in an independent lab. Safe to leave under a patient; similar results as a standard hospital bed sheet.  
(Study available upon request)

Sit-To-Stand



- **One-Size-Fits-Most**
  - Waist size: 30" - 48"
  - Weight capacity: 500 lb (226 kg)
  - Product #: **51961** - same # for VA
  - Sold in cases of 10

One size, fits the average-sized patient, improving compliance and inventory management

Compatible with most loop-style sit-to-stands

Non-slip, padded back panel

Multi-loop adjustment interface

\*Small and large sizes available in reusable option



1-18009-R02