



## **PATIENT SLING COMPATIBILITY STATEMENT**

### Information for Using Wy'East Medical Slings with other Manufacturer's Lift Equipment

Wy'East Medical Corporation has designed and manufactured safe patient handling products since its inception in 1987. Our company is registered with the Food and Drug Administration and is subject to periodic audits by the FDA to verify our compliance with Good Manufacturing Practices. Our slings are sold through national distributors and manufacturers of patient lifting equipment and are routinely used on lift equipment that is not manufactured by Wy'East Medical.

### Lift-Sling Compatibility

Wy'East Medical lift slings are fully compatible with other manufacturers' floor and ceiling lift systems using hook style sling bars. Lifting slings may be used with 2 or 4 point hook style sling bars. Repositioning slings may be used with multiple point hook style sling bars suitable for the size of sling being used. Wy'East Medical floor lifts are fully compatible with other manufacturers' slings utilizing the sling loop attachment method of attaching slings to floor lifts and ceiling lifts. The care provider must follow the responsible vendor's instruction manuals and recommendations with respect to the use, care, cleaning and inspection of slings and patient lifts.

### **COMPATIBLE**

- Loop style Attachment
- Loop attachment on hook style-loop sling compatible sling bar

### **NOT COMPATIBLE**

- Clip style attachment

### Sling Liability

Wy'East Medical stands behind the safe application of our loop slings and accepts responsibility if a product defect results in an error in the practical application of the Wy'East Medical sling.

Wy'East Medical maintains product liability insurance on all the products we manufacture. Wy'East Medical will not be responsible for an error resulting from the improper use or application of the sling and/or lift.