



## Technical Data Sheet

- 51411 Disposable Standard Repositioning Sling, 1000# (1 pkg. - 5 slings)

<b>Device Name</b>	Disposable Standard Repositioning Sling
<b>Device Description</b>	Patient repositioning sling for use with a patient transfer lift and a two point spreader bar to reposition or turn a patient in bed. Intended for single patient use. It is not washable.
<b>Unit Description</b>	1 Disposable Standard Repositioning Sling
<b>Units per Package</b>	5
<b>Package Dimensions (in)</b>	Approximately 17 L x 11 W x 9.5 H (1 pkg. - 5 slings)
<b>Package Weight (lbs)</b>	15
<b>Unit Dimensions (in)</b>	66 L x 42 W
<b>Unit Weight (lbs)</b>	3
<b>Lifting Capacity (lbs)</b>	1000
<b>HIBC Code</b>	M66851411
<b>Materials</b>	Body: Breathable Cotton-Polyester blend. Strapping: Polypropylene Thread: Polyester
<b>Radio Translucent</b>	Yes
<b>Health/Safety</b>	Contains no latex or PVC
<b>Recyclable</b>	No, combustible waste
<b>Packaging</b>	Minimal – all packaging materials are recyclable
<b>FDA Device Class</b>	Class 1
<b>FDA Regulatory Number</b>	880.5510 (component of)
<b>European Device Class</b>	Class 1
<b>European Regulation</b>	Medical Device Directive 93/42/EEC
<b>GMDN Code</b>	40535



Clackamas, Oregon USA ~ [www.wyeastmedical.com](http://www.wyeastmedical.com)  
Toll Free in North America 1.800.255.3126 ~ International +1.503.657.3101