




Reusable Highback Sling

Reusable Slings

Sling Size	51305 Small	51306 Medium	51307 Large	51308 XL
Binding Color	BLUE	MAROON	GREEN	BLACK
Weight Capacity	600 lbs	600 lbs	600 lbs	800 lbs
Patient Weight	45-100 lbs.	100-210 lbs.	210 – 600 lbs.	450 – 800 lbs.
Patient Height	4' - 5'	5' - 6'	6' +	6' +

- * FITS ANY 2-POINT SLING BAR * MACHINE WASHABLE
- * SUPPORTS HEAD * COLOR CODED TO DISTINGUISH SIZES

Patient / Resident Specific Disposable Slings


51170 Disposable Sling 5 pack

51411 Disposable Repositioning Sling 5 pack

- * EASY TO MOVE OR TURN PATIENT * ONE SIZE FITS ALL
- * BREATHABLE— CAN REMAIN UNDER BED'S SHEET
- * STAYS WITH PATIENT FOR ENTIRE HOSPITAL STAY
- * HELPS REDUCE INFECTION AND CONTAMINATION RISK