

SLINGS

All slings fit industry standard
two point sling bars.

Breathable

Rugged

Comfortable

DISPOSABLE HIGH-BACK

Standard Case of 5 (600 lbs.): 51561

Patient Specific Disposable Slings are designed to make it easy to move a patient or resident and can be used with any patient transfer lift. This one-size-fits-all sling can remain under the bed sheet without disrupting patient comfort. The sling can be used multiple times for the same patient or resident. This sling may also be used for lateral transfers. Disposable slings should not be washed.

- Integrated straps for positioning
- Accommodates a wide range of patients




REUSABLE HIGH-BACK

XLarge (800 lbs.): 51308 Medium (600 lbs.): 51306

Large (600 lbs.): 51307 Small (600 lbs.): 51305

The Double Loop High-Back Sling is designed to be used with any loop-style patient transfer lift. The slings stretch slightly to fit the needs of most patients and are clearly labeled to guide head and leg positioning. Color-coded for quick identification: The blue sling is our smallest size, then maroon, green, and black (our largest size, supporting up to 800 pounds). Launder up to 300 times.

- Padded head and leg support
- Color coded for sizing
- Launderable
- Double-loop feature



STAND ASSIST

Small (500 lbs.): 51715 Large (500 lbs.): 51717

Medium (500 lbs.): 51704 Standard Case of 10 (500 lbs.): 51961

Our reusable multi-use belt sling is designed to move patients easily. This sling has a unique 'Sling Management System' which adapts to electronic inventory management systems minimizing the guesswork in sling size selection. This sling incorporates a non-slip padded back panel, multi-loop adjustable interface and is color coded for sizing. Maximum patient weight is 500 pounds.

- Non-slip padded back panel
- Multi-loop adjustment interface
- Color coded for sizing
- Launderable



SLINGS

Multi-loop adjustable head support

Breathable

Rugged

REUSABLE REPOSITIONING

Standard (600 lbs.): 51606

Bariatric (1000 lbs.): 51609

Reusable Repositioning Slings lie flat on the bed and are designed to make it easy to move, or turn a patient, or resident. Composed of breathable material, the sling can remain under the patient without disrupting comfort. These slings are clearly labeled and designed to reduce the risk of infection and contamination. Available in two sizes, these slings can be laundered up to 300 times.

- Launderable



DISPOSABLE REPOSITIONING

Standard Case of 5 (1000 lbs.): 51411

Bariatric (1200 lbs.): 51465

A patient-specific repositioning sling is appropriate in some situations where moving or repositioning patients is required. Composed of breathable material, this sling can remain under the patient without disrupting comfort. These slings are designed to reduce the risk of infection and contamination. They can be used multiple times for the same patient or resident, but they should not be washed or recycled. Available in two sizes.



1.800.255.3126 | fax: 503.657.6901 | P.O. Box 1625 | Clackamas, OR 97015 | www.wyeastmedical.com