



Technical Data Sheet

- 51561 Disposable Double Loop High Back Patient Lifting Sling (1 pkg. - 5 slings)

Device Name	Disposable Double Loop High Back Patient Lifting Sling
Device Description	Patient transfer sling for use with a patient transfer lift and a two point spreader bar to move a patient between bed, chair, commode, floor and other locations. Intended for single patient use. It is not washable.
Unit Description	1 Disposable Double Loop High Back Patient Lifting Sling
Units per Package	5
Package Dimensions (in)	Approximately 14 L x 12 W x 7.5 H (1 pkg. - 5 slings)
Package Weight (lbs)	9
Unit Dimensions (in)	77 L x 44 W
Unit Weight (lbs)	1.8
Lifting Capacity (lbs)	600
HIBC Code	M66851561
Materials	Body: Breathable Cotton - Polyester blend. Color: White Strapping: Polypropylene Thread: Polyester
Radio Translucent	Yes
Health/Safety	Contains no latex or PVC
Recyclable	No, combustible waste
Packaging	Minimal – all packaging materials are recyclable
FDA Device Class	Class 1
FDA Regulatory Number	880.5510 (component of)
European Device Class	Class 1
European Regulation	Medical Device Directive 93/42/EEC
GMDN Code	40535



Clackamas, Oregon USA ~ www.wyeastmed.com

Toll Free in North America 1.800.255.3126 ~ International +1.503.657.3101